

Purified anti-human IL-12/IL-23 p40

Mouse IgG1, kappa
C11.5



BioLegend
From Revvity

500 µg
Artikel-Nr.: B501802 | **BioLegend** | Hersteller-Nr.: 501802

207,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: C11.5
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human IL-12/IL-23 p40

Clone: C11.5

Isotype: Mouse IgG1, kappa

Other_Names: Interleukin-12 p40, Interleukin-23 p40, Cytotoxic lymphocyte maturation factor (CLMF), Natural killer cell stimulatory factor (NKSF), CTL maturation factor (TcMF), T-cell stimulating factor (TSF)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1023>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>