

Purified anti-human IL-13

Rat IgG1, kappa
JES10-5A2



BioLegend
From Revvity

50 µg
Artikel-Nr.: B501901 | [BioLegend](#) | Hersteller-Nr.: 501901

**85,00 €*
zzgl. Versandkosten**

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JES10-5A2
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-human IL-13

Clone: JES10-5A2
Isotype: Rat IgG1, kappa
Other_Names: Interleukin-13, NC30
Size: 50 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1338>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>