

Purified anti-human G-CSF

Rat IgG1, kappa
BVD13-3A5



BioLegend
From Revvity

500 µg

Artikel-Nr.: B502102 | [BioLegend](#) | Hersteller-Nr.: 502102

188,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: BVD13-3A5
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-human G-CSF

Clone: BVD13-3A5

Isotype: Rat IgG1, kappa

Other_Names: Granulocyte-Colony Stimulating Factor, CSF-beta, Pluripoietin, Granulocytic neutrophil colony stimulating activity (G-CSA), Macrophage-granulocyte inducer-2 (MGI-2)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=914>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>