

Purified anti-human IFN-gamma

Mouse IgG1, kappa
4S.B3



BioLegend
From Revvity

500 µg
Artikel-Nr.: B502502 | **BioLegend** | Hersteller-Nr.: 502502

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 4S.B3
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human IFN-gamma

Clone: 4S.B3

Isotype: Mouse IgG1, kappa

Other_Names: Interferon-gamma, Immune interferon, Type II interferon, T cell interferon, Macrophage-activating factor (MAF), IFN-g, IFN-gamma

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1015>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>