

Alexa Fluor® 647 anti-human IL-4

Rat IgG1, kappa
MP4-25D2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B500818 | [BioLegend](#) | Hersteller-Nr.: 500818

212,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: MP4-25D2
Format: A647
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human IL-4

Clone: MP4-25D2

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-4, Ia inducing factor (IaIF), B-cell stimulating factor-1 (BSF-1), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 100 tests

Reactivity: Human

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2745>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>