

## Purified anti-human IL-5

Rat IgG2a, kappa  
JES1-39D10



**BioLegend**  
From Revvity

500 µg  
Artikel-Nr.: B500902 | **BioLegend** | Hersteller-Nr.: 500902

**210,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** JES1-39D10  
**Format:** PURE  
**Isotype:** Rat IgG2a, kappa  
**Species of Origin:** Rat  
**Verpackung:** 500 µg

Purified anti-human IL-5

Clone: JES1-39D10

Isotype: Rat IgG2a, kappa

Other\_Names: Interleukin-5, Eosinophil differentiation factor (EDF), Eosinophil colony stimulating factor (Eo-CSF), B cell growth factor-2 (BCGF-2), B cell differentiation factor for IgM (BCDF-m), IgA enhancing factor (IgA-EF), T cell replacing factor-1 (TRF-1)

Size: 500 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=935>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>