

PE anti-human IL-5

Rat IgG2a, kappa
JES1-39D10



BioLegend
From Revvity

100 tests

Artikel-Nr.: B500904 | [BioLegend](#) | Hersteller-Nr.: 500904

196,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JES1-39D10
Format: PE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 Tests

PE anti-human IL-5

Clone: JES1-39D10

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-5, Eosinophil differentiation factor (EDF), Eosinophil colony stimulating factor (Eo-CSF), B cell growth factor-2 (BCGF-2), B cell differentiation factor for IgM (BCDF-m), IgA enhancing factor (IgA-EF), T cell replacing factor-1 (TRF-1)

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=934>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>