

Purified anti-human IL-6

Rat IgG1, kappa
MQ2-13A5



BioLegend
From Revvity

50 µg
Artikel-Nr.: B501101 | [BioLegend](#) | Hersteller-Nr.: 501101

67,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MQ2-13A5
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-human IL-6

Clone: MQ2-13A5

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-6, B cell stimulating factor-2 (BSF-2), Cytotoxic T cell differentiation factor (CDF), Hepatocyte stimulating factor (HSF), Hybridoma/plasmacytoma growth factor (HPGF)

Size: 50 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=983>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>