

Alexa Fluor® 647 anti-human IL-10

Rat IgG1, kappa
JES3-9D7



BioLegend
From Revvity

100 tests

Artikel-Nr.: B501412 | [BioLegend](#) | Hersteller-Nr.: 501412

212,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JES3-9D7
Format: A647
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human IL-10

Clone: JES3-9D7

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-10, B cell derived T cell growth factor (B-TCGF), Cytokine synthesis inhibitory factor (CSIF), T-cell growth inhibitory factor (TGIF)

Size: 100 tests

Reactivity: Human

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2747>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>