

FITC anti-human IL-4

Rat IgG1, kappa
MP4-25D2



BioLegend
From Revvity

50 µg
Artikel-Nr.: B500807 | **BioLegend** | Hersteller-Nr.: 500807

100,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MP4-25D2
Format: FITC
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 50 µg

FITC anti-human IL-4

Clone: MP4-25D2

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-4, Ia inducing factor (IaIF), B-cell stimulating factor-1 (BSF-1), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 50 µg

Reactivity: Human

Format: FITC

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=967>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>