

## Purified anti-human IL-1alpha

Mouse IgG1, kappa  
364-3B3-14



**BioLegend**  
From Revvity

500 µg

Artikel-Nr.: B500104 | **BioLegend** | Hersteller-Nr.: 500104

**216,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** 364-3B3-14  
**Format:** PURE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 500 µg

Purified anti-human IL-1alpha

Clone: 364-3B3-14

Isotype: Mouse IgG1, kappa

Other\_Names: Interleukin-1alpha, Lymphocyte-activating factor (LAF), Endogenous pyrogen (EP), Leukocyte endogenous mediator (LEM), Mononuclear cell factor (MCF/MNCF)

Size: 500 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1325>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>