

PE anti-human IL-1alpha

Mouse IgG1, kappa
364-3B3-14



BioLegend
From Revvity

100 tests

Artikel-Nr.: B500106 | **BioLegend** | Hersteller-Nr.: 500106

244,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 364-3B3-14
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human IL-1alpha

Clone: 364-3B3-14

Isotype: Mouse IgG1, kappa

Other_Names: Interleukin-1alpha, Lymphocyte-activating factor (LAF), Endogenous pyrogen (EP), Leukocyte endogenous mediator (LEM), Mononuclear cell factor (MCF/MNCF)

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1324>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>