

Purified anti-human IL-2

Rat IgG2a, kappa
MQ1-17H12



BioLegend
From Revvity

500 µg
Artikel-Nr.: B500302 | **BioLegend** | Hersteller-Nr.: 500302

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MQ1-17H12
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-human IL-2

Clone: MQ1-17H12

Isotype: Rat IgG2a, kappa

Other_Names: Interleukin-2, T cell growth factor (TCGF), Killer cell helper factor (KHF), Macrophage-activating factor for cytotoxicity I (MAF-C I)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1352>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>