

APC anti-human IL-2

Rat IgG2a, kappa
MQ1-17H12



BioLegend
From Revvity

50 µg
Artikel-Nr.: B500311 | [BioLegend](#) | Hersteller-Nr.: 500311

133,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: MQ1-17H12
Format: APC
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 50 µg

APC anti-human IL-2

Clone: MQ1-17H12
Isotype: Rat IgG2a, kappa
Other_Names: Interleukin-2, T cell growth factor (TCGF), Killer cell helper factor (KHF), Macrophage-activating factor for cytotoxicity I (MAF-C I)
Size: 50 µg

Reactivity: Human
Format: APC
App_Abbrev: FC, ICFC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1348>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>