

PE anti-human IL-3

Rat IgG1, kappa
BVD3-1F9



BioLegend
From Revvity

100 tests

Artikel-Nr.: B500606 | **BioLegend** | Hersteller-Nr.: 500606

244,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: BVD3-1F9
Format: PE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 Tests

PE anti-human IL-3

Clone: BVD3-1F9

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-3, Burst promoting activity, Eosinophil colony stimulating factor (Eo-CSF), Hematopoietic cell growth factor (HCGF), Mast (MGF/MCGF), Multi-colony stimulating (Multi-CSF), P cell stimulating activity (PCSA), Thy1 inducing factor

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=921>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>