

Purified anti-mouse Ig light chain lambda

Rat IgG2a, kappa
RML-42



BioLegend
From Revvity

500 µg
Artikel-Nr.: B407302 | **BioLegend** | Hersteller-Nr.: 407302

137,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: RML-42
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse Ig light chain lambda

Clone: RML-42
Isotype: Rat IgG2a, kappa
Other_Names: Immunoglobulin light chain lambda
Size: 500 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2631>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>