

Purified anti-rat IgG2c

Mouse IgG1, kappa
MRG2c-67



BioLegend
From Revvity

500 µg

Artikel-Nr.: B407602 | **BioLegend** | Hersteller-Nr.: 407602

137,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MRG2c-67
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-rat IgG2c

Clone: MRG2c-67
Isotype: Mouse IgG1, kappa
Other_Names: Immunoglobulin G2c
Size: 500 µg

Reactivity: Rat
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2637>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>