

Purified anti-mouse IgMb

Mouse (BALB/c) IgG1, kappa
AF6-78



BioLegend
From Revvity

500 µg

Artikel-Nr.: B406202 | **BioLegend** | Hersteller-Nr.: 406202

207,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA, FC
Clone: AF6-78
Format: PURE
Isotype: Mouse (BALB/c) IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-mouse IgMb

Clone: AF6-78
Isotype: Mouse (BALB/c) IgG1, kappa
Other_Names: Igh-6b
Size: 500 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ELISA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1746>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>