

Purified anti-mouse IgG1

Rat IgG
RMG1-1



BioLegend
From Revvity

50 µg
Artikel-Nr.: B406601 | [BioLegend](#) | Hersteller-Nr.: 406601

57,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: RMG1-1
Format: PURE
Isotype: Rat IgG
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse IgG1

Clone: RMG1-1
Isotype: Rat IgG
Other_Names: Immunoglobulin G1
Size: 50 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2516>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>