

Purified anti-mouse IgA

Rat IgG1, kappa
RMA-1



BioLegend
From Revvity

500 µg

Artikel-Nr.: B407002 | **BioLegend** | Hersteller-Nr.: 407002

137,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: RMA-1
Format: PURE
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 500 µg

Purified anti-mouse IgA

Clone: RMA-1
Isotype: Rat IgG1, kappa
Other_Names: Immunoglobulin A
Size: 500 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2519>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>