

Purified Goat anti-mouse IgG (minimal x-reacti...

Goat Polyclonal IgG
Poly4053



BioLegend
From Revvity

500 µg
Artikel-Nr.: B405301 | **BioLegend** | Hersteller-Nr.: 405301

146,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: Poly4053
Format: PURE
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 500 µg

Purified Goat anti-mouse IgG (minimal x-reactivity)

Clone: Poly4053
Isotype: Goat Polyclonal IgG
Other_Names:
Size: 500 µg

Reactivity: Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1419>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>