

Biotin Goat anti-mouse IgG (minimal x-reactivi...

Goat Polyclonal IgG
Poly4053



BioLegend
From Revvity

500 µg
Artikel-Nr.: B405303 | **BioLegend** | Hersteller-Nr.: 405303

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: Poly4053
Format: BIOTIN
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 500 µg

Biotin Goat anti-mouse IgG (minimal x-reactivity)

Clone: Poly4053
Isotype: Goat Polyclonal IgG
Other_Names:
Size: 500 µg

Reactivity: Mouse
Format: BIOTIN
App_Abbrev: FC, ICFC
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1387>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>