

PE/Cyanine5 Rat IgG2b, kappa Isotype Ctrl

Rat IgG2b, kappa
RTK4530



BioLegend
From Revvity

25 µg
Artikel-Nr.: B400609 | **BioLegend** | Hersteller-Nr.: 400609

76,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: RTK4530
Format: PE_CYANINE5
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 25 µg

PE/Cyanine5 Rat IgG2b, kappa Isotype Ctrl

Clone: RTK4530
Isotype: Rat IgG2b, kappa
Other_Names:
Size: 25 µg

Reactivity:
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1857>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>