

Purified Mouse IgG2b, kappa Isotype Ctrl

Mouse IgG2b, kappa
MPC-11



BioLegend
From Revvity

50 µg
Artikel-Nr.: B400301 | **BioLegend** | Hersteller-Nr.: 400301

24,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MPC-11
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 50 µg

Purified Mouse IgG2b, kappa Isotype Ctrl

Clone: MPC-11
Isotype: Mouse IgG2b, kappa
Other_Names:
Size: 50 µg

Reactivity:
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1416>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>