

PE/Cyanine5 Mouse IgG2b, kappa Isotype Ctrl

Mouse IgG2b, kappa
MPC-11



BioLegend
From Revvity

25 tests

Artikel-Nr.: B400317 | [BioLegend](#) | Hersteller-Nr.: 400317

85,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: MPC-11
Format: PE_CYANINE5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE/Cyanine5 Mouse IgG2b, kappa Isotype Ctrl

Clone: MPC-11
Isotype: Mouse IgG2b, kappa
Other_Names:
Size: 25 tests

Reactivity:
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1415>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>