

PE/Cyanine5 Mouse IgG2a, kappa Isotype Ctrl

Mouse IgG2a, kappa
MOPC-173



BioLegend
From Revvity

100 tests

Artikel-Nr.: B400218 | BioLegend | Hersteller-Nr.: 400218

132,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MOPC-173
Format: PE_CYANINE5
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 Mouse IgG2a, kappa Isotype Ctrl

Clone: MOPC-173
Isotype: Mouse IgG2a, kappa
Other_Names:
Size: 100 tests

Reactivity:
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1402>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>