

Purified anti-human CD337 (NKp30)

Mouse IgG1, kappa
P30-15



BioLegend
From Revvity

100 µg

Artikel-Nr.: B325202 | [BioLegend](#) | Hersteller-Nr.: 325202

104,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Block, ELISA, FC
Clone: P30-15
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD337 (NKp30)

Clone: P30-15

Isotype: Mouse IgG1, kappa

Other_Names: NKp30, NCR3, Activating NK receptor NKp30, natural cytotoxicity triggering receptor 3

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ELISA, Block

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3852>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>