

PE anti-human CD337 (NKp30)

Mouse IgG1, kappa
P30-15



BioLegend
From Revvity

25 tests

Artikel-Nr.: B325207 | **BioLegend** | Hersteller-Nr.: 325207

108,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: P30-15
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human CD337 (NKp30)

Clone: P30-15

Isotype: Mouse IgG1, kappa

Other_Names: NKp30, NCR3, Activating NK receptor NKp30, natural cytotoxicity triggering receptor 3

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3855>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>