

APC anti-human CD337 (NKp30)

Mouse IgG1, kappa
P30-15



BioLegend
From Revvity

100 tests

Artikel-Nr.: B325210 | BioLegend | Hersteller-Nr.: 325210

252,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: P30-15
Format: APC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

APC anti-human CD337 (NKp30)

Clone: P30-15

Isotype: Mouse IgG1, kappa

Other_Names: NKp30, NCR3, Activating NK receptor NKp30, natural cytotoxicity triggering receptor 3

Size: 100 tests

Reactivity: Human

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3856>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>