

## Alexa Fluor® 647 anti-human CD337 (NKp30)

Mouse IgG1, kappa  
P30-15



**BioLegend**  
From Revvity

100 tests

Artikel-Nr.: B325212 | BioLegend | Hersteller-Nr.: 325212

**235,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** P30-15  
**Format:** A647  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 Tests

Alexa Fluor® 647 anti-human CD337 (NKp30)

Clone: P30-15

Isotype: Mouse IgG1, kappa

Other\_Names: NKp30, NCR3, Activating NK receptor NKp30, natural cytotoxicity triggering receptor 3

Size: 100 tests

Reactivity: Human

Format: A647

App\_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3857>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>