

Alexa Fluor® 647 anti-human CD326 (EpCAM)

Mouse IgG2b, kappa
9C4



BioLegend
From Revvity

100 tests

Artikel-Nr.: B324212 | [BioLegend](#) | Hersteller-Nr.: 324212

235,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 9C4
Format: A647
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human CD326 (EpCAM)

Clone: 9C4

Isotype: Mouse IgG2b, kappa

Other_Names: Ep-CAM, tumor associated calcium signal transducer 1, epithelial cell surface antigen, epithelial glycoprotein 2, EGP2, adenocarcinoma associated antigen, TROP1

Size: 100 tests

Reactivity: Human

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3760>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>