

Purified anti-human CD334 (FGFR4)

Mouse IgG1, kappa
4FR6D3



BioLegend
From Revvity

100 µg

Artikel-Nr.: B324302 | BioLegend | Hersteller-Nr.: 324302

165,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 4FR6D3
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD334 (FGFR4)

Clone: 4FR6D3

Isotype: Mouse IgG1, kappa

Other_Names: FGFR4, fibroblast growth factor receptor 4, TKF, FGF R4

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3761>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>