

PE anti-human CD334 (FGFR4)

Mouse IgG1, kappa
4FR6D3



BioLegend
From Revvity

100 tests

Artikel-Nr.: B324306 | BioLegend | Hersteller-Nr.: 324306

234,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 4FR6D3
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD334 (FGFR4)

Clone: 4FR6D3
Isotype: Mouse IgG1, kappa
Other_Names: FGFR4, fibroblast growth factor receptor 4, TKF, FGF R4
Size: 100 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3763>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>