

Purified anti-human CD340 (erbB2/HER-2)

Mouse IgG1, kappa
24D2



BioLegend
From Revvity

25 µg
Artikel-Nr.: B324401 | **BioLegend** | Hersteller-Nr.: 324401

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 24D2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD340 (erbB2/HER-2)

Clone: 24D2

Isotype: Mouse IgG1, kappa

Other_Names: HER2/erbB2, tyrosine kinase cell surface receptor HER2, oncogene ERBB2, oncogene NGL, neuroblastoma- or glioblastoma-derived, HER-2

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3764>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>