

PE anti-human erbB3/HER-3

Mouse IgG2a, kappa
1B4C3



BioLegend
From Revvity

25 tests

Artikel-Nr.: B324705 | **BioLegend** | Hersteller-Nr.: 324705

108,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 1B4C3
Format: PE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human erbB3/HER-3

Clone: 1B4C3
Isotype: Mouse IgG2a, kappa
Other_Names: erbB-3, erbB3, HER-3, HER3
Size: 25 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3745>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>