

APC anti-human CD336 (NKp44)

Mouse IgG1, kappa
P44-8



BioLegend
From Revvity

25 tests

Artikel-Nr.: B325109 | BioLegend | Hersteller-Nr.: 325109

120,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: P44-8
Format: APC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

APC anti-human CD336 (NKp44)

Clone: P44-8

Isotype: Mouse IgG1, kappa

Other_Names: NKp44, NCR2, Activating NK receptor NKp44, natural cytotoxicity triggering receptor 2. lymphocyte antigen 95 homolog

Size: 25 tests

Reactivity: Human

Format: APC

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3850>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>