

Alexa Fluor® 647 anti-human CD336 (NKp44)

Mouse IgG1, kappa
P44-8



BioLegend
From Revvity

100 tests

Artikel-Nr.: B325112 | BioLegend | Hersteller-Nr.: 325112

238,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: P44-8
Format: A647
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human CD336 (NKp44)

Clone: P44-8

Isotype: Mouse IgG1, kappa

Other_Names: NKp44, NCR2, Activating NK receptor NKp44, natural cytotoxicity triggering receptor 2. lymphocyte antigen 95 homolog

Size: 100 tests

Reactivity: Human

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3851>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>