

Purified anti-human CD140a (PDGFRalpha)

Mouse IgG1, kappa
16A1



BioLegend
From Revvity

100 µg
Artikel-Nr.: B323502 | **BioLegend** | Hersteller-Nr.: 323502

95,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 16A1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD140a (PDGFRalpha)

Clone: 16A1

Isotype: Mouse IgG1, kappa

Other_Names: Platelet-derived growth factor receptor, alpha polypeptide, PDGFR2, PDGFRalpha, PDGFRa, PDGF receptor alpha

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3725>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>