

PE anti-human CD140b (PDGFRbeta)

Mouse IgG1, kappa
18A2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B323606 | [BioLegend](#) | Hersteller-Nr.: 323606

215,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 18A2
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD140b (PDGFRbeta)

Clone: 18A2

Isotype: Mouse IgG1, kappa

Other_Names: Platelet-derived growth factor receptor, beta polypeptide, PDGFR1, PDGFRbeta, PDGFRb, PDGF receptor beta

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3730>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>