

PE anti-human CD172b (SIRPbeta)

Mouse IgG1, kappa
B4B6



BioLegend
From Revvity

100 tests

Artikel-Nr.: B323906 | **BioLegend** | Hersteller-Nr.: 323906

210,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: B4B6
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD172b (SIRPbeta)

Clone: B4B6
Isotype: Mouse IgG1, kappa
Other_Names: SIRPb, CD172 antigen-like family member B
Size: 100 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3748>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>