

FITC anti-human CD318 (CDCP1)

Mouse IgG2b, kappa
CUB1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B324004 | BioLegend | Hersteller-Nr.: 324004

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: CUB1
Format: FITC
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

FITC anti-human CD318 (CDCP1)

Clone: CUB1
Isotype: Mouse IgG2b, kappa
Other_Names: CDCP1, CUB domain containing protein 1, SIMA135
Size: 100 tests

Reactivity: Human
Format: FITC
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3732>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>