

PE anti-human TSLPR (TSLP-R)

Mouse IgG1, kappa
1B4



BioLegend
From Revvity

25 tests

Artikel-Nr.: B322805 | **BioLegend** | Hersteller-Nr.: 322805

113,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 1B4
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human TSLPR (TSLP-R)

Clone: 1B4

Isotype: Mouse IgG1, kappa

Other_Names: Thymic stromal lymphopoietin protein receptor, cytokine receptor-like 2, CRL2, CRLF2, IL-XR, TSLP-R

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3692>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>