

Biotin anti-human CD47

Mouse IgG1, kappa
CC2C6



BioLegend
From Revvity

100 µg

Artikel-Nr.: B323104 | BioLegend | Hersteller-Nr.: 323104

169,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: CC2C6
Format: BIOTIN
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Biotin anti-human CD47

Clone: CC2C6

Isotype: Mouse IgG1, kappa

Other_Names: Rh-associated protein, gp42, integrin-associated protein, IAP, neurophilin

Size: 100 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3706>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>