

Purified anti-human CD117 (c-kit)

Mouse IgG1, kappa
A3C6E2



BioLegend
From Revvity

100 µg

Artikel-Nr.: B323402 | **BioLegend** | Hersteller-Nr.: 323402

141,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: A3C6E2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD117 (c-kit)

Clone: A3C6E2

Isotype: Mouse IgG1, kappa

Other_Names: Stem cell factor receptor, c-kit, mast cell growth factor receptor, steel factor receptor

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3715>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>