

Biotin anti-human CD117 (c-kit)

Mouse IgG1, kappa
A3C6E2



BioLegend
From Revvity

100 µg

Artikel-Nr.: B323406 | [BioLegend](#) | Hersteller-Nr.: 323406

165,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: A3C6E2
Format: BIOTIN
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Biotin anti-human CD117 (c-kit)

Clone: A3C6E2

Isotype: Mouse IgG1, kappa

Other_Names: Stem cell factor receptor, c-kit, mast cell growth factor receptor, steel factor receptor

Size: 100 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3717>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>