

## PE anti-human CD117 (c-kit)

Mouse IgG1, kappa  
A3C6E2



**BioLegend**  
From Revvity

100 tests

Artikel-Nr.: B323408 | BioLegend | Hersteller-Nr.: 323408

**215,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** A3C6E2  
**Format:** PE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 Tests

PE anti-human CD117 (c-kit)

Clone: A3C6E2

Isotype: Mouse IgG1, kappa

Other\_Names: Stem cell factor receptor, c-kit, mast cell growth factor receptor, steel factor receptor

Size: 100 tests

Reactivity: Human

Format: PE

App\_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3718>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>