

Purified anti-human Lymphotoxin beta Receptor ...

Mouse IgG2b, kappa
31G4D8



BioLegend
From Revvity

100 µg

Artikel-Nr.: B322002 | [BioLegend](#) | Hersteller-Nr.: 322002

113,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 31G4D8
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human Lymphotoxin beta Receptor (LTbetaR)

Clone: 31G4D8

Isotype: Mouse IgG2b, kappa

Other_Names: Lymphotoxin beta receptor, Tumor necrosis factor receptor superfamily member 3

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3282>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>