

PE/Cyanine5 anti-human CD206 (MMR)

Mouse IgG1, kappa
15-2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B321108 | [BioLegend](#) | Hersteller-Nr.: 321108

244,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 15-2
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD206 (MMR)

Clone: 15-2

Isotype: Mouse IgG1, kappa

Other_Names: MMR (macrophage mannose receptor), MR (mannose receptor), CD206, MRC1

Size: 100 tests

Reactivity: Human

Format: PE_CYANINE5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2995>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>