

PE anti-human/mouse integrin beta7

Rat IgG2a, kappa
FIB504



BioLegend
From Revvity

25 tests

Artikel-Nr.: B321203 | **BioLegend** | Hersteller-Nr.: 321203

98,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: FIB504
Format: PE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 25 tests

PE anti-human/mouse integrin beta7

Clone: FIB504
Isotype: Rat IgG2a, kappa
Other_Names: beta7 Integrin, integrin betap, ITGB7
Size: 25 tests

Reactivity: Human, Mouse, Cynomolgus, Rhesus
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2952>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>