

Purified anti-human/mouse CLA

Rat IgM, kappa
HECA-452



BioLegend
From Revvity

100 µg

Artikel-Nr.: B321302 | **BioLegend** | Hersteller-Nr.: 321302

96,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: HECA-452
Format: PURE
Isotype: Rat IgM, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human/mouse CLA

Clone: HECA-452
Isotype: Rat IgM, kappa
Other_Names: Cutaneous Lymphocyte-associated Antigen (CLA)
Size: 100 µg

Reactivity: Human, Mouse
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2962>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>